

COMPARISON WITH EXISTING SOLUTIONS

NOBLE DOME:

- IS the glove box
- No sample size limit
- No stage travel limit
- No mechanic or electronic parts to break or repair
- Air-free TEM rod easily fits inside, streamlining FIB to TEM sample transfers

OTHER PRODUCTS:

- Assumes users have their own glove box (\$10-100K cost)
- Limits sample size
- Limits stage travel
- Overly intricate: self-service difficult & can malfunction inside instrument
- Complex workflow for TEM preparation



THE ONLY GLOVE BOX LOAD-LOCK FOR ELECTRON MICROSCOPY

An innovative solution for loading air-sensitive materials, eliminating the need for an intermediary device to transfer samples between instruments



Transparent dome lets users stand directly over samples for easier manipulation



THE NOBLE DOME WAS INVENTED & PROTOTYPED AT THE CENTER FOR ADVANCED MATERIALS CHARACTERIZATION OREGON, UNIVERSITY OF OREGON WITH A PATENT PENDING



CONTACT INFO: VALERIE BROGDEN
valerie@nobledome.net
nobledome.net

**AIR-FREE
TRANSFER SYSTEM**


UNIQUE

NOBLE DOME

TRADITIONAL LOADING PROCEDURES EXPOSE SAMPLES TO ATMOSPHERE DURING TRANSFER, THE NOBLE DOME PRESERVES THE INTEGRITY OF AIR-SENSITIVE MATERIALS BY TRANSFERRING THEM IN AND OUT OF AN INSTRUMENT IN AN OXYGEN-FREE ENVIRONMENT

SIMPLIFY YOUR WORKFLOW:

CONVENTIONAL METHOD: COMPLICATED WITH MULTIPLE DEVICES & STANDALONE GLOVE BOX




NOBLE DOME APPROACH:
LOAD SAMPLES DIRECTLY
ONTO TEM HOLDER FROM
THE SEM/FIB

LITHIUM METAL OXIDATION EXPERIMENT:

The same location on a lithium metal sample was imaged & analyzed with EDS in its pristine state, after 30 minutes in an argon environment in Noble Dome, & after 30 minutes in atmosphere

Very little change seen in surface topography & oxygen counts between sample in its pristine state & after being in Noble Dome
Notable increase in surface oxidation & topography observed on sample after sitting in atmosphere






PRISTINE SAMPLE

30 MINUTES IN NOBLE DOME

30 MINUTES IN AIR

SEM IMAGES



OXYGEN EDS MAPS

each yellow pixel represents one count

